



Optec s.c.

ul. Alberta 6E, 16-001 Księżyno k. Białegostoku

tel. (85) 743 74 17, (85) 878 16 10

e-mail: medycyna@optec.pl



Ginekologia



KAMERA UC-100 FULL HD USB 3.0 + LED

www.medycyna.sklep.pl



Ginekologia



GQ-6 Hysteroscopy set

F0107 0° Hysteroscopy
Ø6,6x227mm
F0108 0° Hysteroscopy
Ø5,3x224mm

F0101.1 22° Hysteroscopy
Ø6,6x175mm
F0101.3 22° Hysteroscopy
Ø6,6x200mm
F0106 22° Hysteroscopy
Ø5,3x170mm



1. Connecting to TIANSONG light guide cable OLYMPUS
Connecting to OLYMPUS light guide cable STORZ
2. Connecting to STORZ light guide cable WOLF
3. Connecting to WOLF light guide cable ACMI, STRYKER
Connecting to ACMI, STRYKER light guide cable



Flexible operating forceps Ø2,4x395mm



F2101 Biopsy forceps




F2102 Scissor



F2107 Loop catching forceps



F2108 Scissor



F2103 Grasping forceps with 4 claws



Flexible operating forceps Ø1,8x395mm, 5,3mm



F1110 Biopsy forceps



F1111 Foreign body forceps



F1112 Serrated forceps




F1113 Scissors



F1114 Grasping forceps with 4 claws



F2104 Bush right Ø2,5x360mm
F2104.1 Bush flexible Ø1,8x320mm



F4003 Irrigation pump 0.1.1
Flow rate: 0.1L/min



F4001.1 Infusion device (Disposable)



T5010 Light guide cable Ø4,5x2000mm



A2145 Sealing cap



Ginekologia

GQ-3 Hysteroscopy set




A0012 30° Endoscope
Ø4x302mm




F1115 24Fr Operating sheath

GQ-13 Hysteroscopy set




A0012 30° Endoscope
Ø4x302mm



F1116.3 Outer sheath 30° 21r


F1117.3 Manipulator

GQ-18 Hysteroscopy set




F 0012 30° Endoscope 16,5Fr x 201mm

GQ-3, 13, 18 Hysteroscopy set interchangeable parts



F2110 Electric coagulator, Blunt head Ø1,8x410mm



Endoscopy forceps



A2101.1 Biopsy forceps, rigid
7Fr x 395mm



A2102.1 Foreign body forceps, rigid
7Fr x 395mm



A2104.1 Scissors, rigid
7Fr x 395mm



A2107.3 Loop catching forceps, rigid
7Fr x 395mm



M1002.1 High frequency cord, Monopolar



T5010 Light guide cable Ø4,5x2000mm



T1001 Adaptor



T1009.1 Luer-lock adaptor



Ginekologia



GQ-8 Hysteroscopy set (Examine)



T1009.1 Luer-lock adaptor



T1001 Adaptor



T5010 Light guide cable $\varnothing 4,5 \times 2000 \text{ mm}$



A0012 30° Endoscope $\varnothing 4 \times 302 \text{ mm}$



F1101.1 Outer sheath 17,5Fr



F1102.1 Inner sheath 14Fr

GQ-10 Hysteroscopy set (Examine)




T1009.1 Luer-lock adaptor



T1001 Adaptor



T5010 Light guide cable $\varnothing 4,5 \times 2000 \text{ mm}$



A0022 Endoscope $\varnothing 2,9 \times 302 \text{ mm}$




F1104 Outer sheath 15Fr



F1105 Inner sheath 11,7Fr


GQ-12 Hysteroscopy set





Ginekologia


GQ-15 Hysteroscopy set



GQ-16 Hysteroscopy set



GQ-12, 15, 16 Hysteroscopy set interchangeable parts





Ginekologia




DQ-5 Resectoscopy set



A0011 12° Endoscope
Ø4x302mm

1. OLYMPUS
STORZ
2. WOLF
3. ACMI, STRYKER

1. Connecting to TIANSONG light guide cable OLYMPUS
Connecting to OLYMPUS light guide cable STORZ
2. Connecting to STORZ light guide cable WOLF
Connecting to WOLF light guide cable ACMI, STRYKER
3. Connecting to ACMI, STRYKER light guide cable



A0011 12° Endoscope Ø4x302mm



A3162 Working element (Passive)



A3163.2 Inner sheath with obturator 22Fr



A3165.2 Outer sheath 25,5Fr



Electrode



A3111 Loop



A3122 Roller



A3114 Straight loop



T1009.1 Luer-lock adaptor



T1001 Adaptor



M0015 H.F. Generator (Multifunction)



M1001 High frequency cord (Monopolar)



T5010 Light guide cable Ø4,5x2000mm



Ginekologia

Electric hystera-cutter



F5002 Control unit

F3002.2 Hand piece



Abdominal cutter
F2012.1 $\varnothing 10 \times 260 \text{ mm}$
F2012.2 $\varnothing 15 \times 260 \text{ mm}$
F2012.7 $\varnothing 18 \times 260 \text{ mm}$



Trocarn cannula
F2023 $\varnothing 15 \text{ mm}$
F2026 $\varnothing 18 \text{ mm}$



Vagina cutter
F2012.3 $\varnothing 18 \times 260 \text{ mm}$
F2012.6 $\varnothing 15 \times 260 \text{ mm}$



Dilator
F2024 $\varnothing 10 \text{ mm} / \varnothing 15 \text{ mm}$
F2028 $\varnothing 10 \text{ mm} / \varnothing 18 \text{ mm}$



F2011 Uterine moving bar $\varnothing 5 \times 450 \text{ mm}$



Dilator
F3025 $\varnothing 15 / \varnothing 10 \text{ mm}$
F3025.1 $\varnothing 15 / \varnothing 5 \text{ mm}$
F3027 $\varnothing 18 / \varnothing 10 \text{ mm}$
F3031 $\varnothing 10 / \varnothing 5 \text{ mm}$



F2016 Guiding bar 10mm



Hysteromyoma punch
F2017 10mm **F2018** 5mm




F2010 Hysterometer $\varnothing 5 \times 550 \text{ mm}$




Ginekologia



Electric hystera-cutter



F2009 Uterine manipulator



F2005 Salpinx Grasping forceps
Ø5x400mm

F2008.2 Canine Biopsy forceps
Ø5x400mm

F2006.1 Uterine grasping forceps
Ø5x400mm




F2006 Uterine grasping forceps
Ø10x400mm

F2006.3 Uterine grasping forceps
Ø5x400mm



Cervical forceps

F2007 Single claw




Cervical forceps

F2007.1 Double claw



Cervical forceps


F2007.3 Double claw, curved



Simple uterine manipulator 250mm

F2009.2 Double claw

Gynecology instruments



Vaginal tractor

F2029.2 27x70mm
F2029.3 30x80mm
F2029.4 36x90mm
F2029.6 36x110mm
F2029.7 30x110mm
F2029.8 27x110mm
F2029.9 21x110mm
F2029.10 15x110mm

Vaginal tractor

F2040.1 130x20mm
F2040.2 130x25mm
F2040.3 130x30mm
F2040.4 130x35mm
F2040.5 130x40mm
F2040.6 160x30mm
F2040.7 160x40mm
F2040.8 160x45mm

Vaginal tractor

F2041.1 45x115mm
F2041.2 45x85mm
F2041.3 45x55mm



Ginekologia

Gynecology instruments




F2030 320mm Hysteromyoma dissector




F2035 250mm Cervical plate




F2036 250mm Vaginal plate



F2032 Ligamenta hook forceps



F2034 Myoma grasping forceps



F2043 Tissue forceps, curved, loop catching



F2044 Grasping forceps



Vaginal speculum
F2042.1 35x100mm
F2042.2 35x85mm
F2042.3 32x75mm
F2042.5 17x75mm



F2004 Oviduct ligature forceps
ø5,8x350mm







Ginekologia



Gynecology instruments





Ginekologia

Gynecology instruments



F2033.2 Cup uterine manipulator, Autoclave



Cup 1



Cup 2



Cup 3



Holder



Holder



Head



Head



F2033.3 Uterine manipulator, Distraction, curved



Ginekologia



Endoscopy camera system

4K Endoscope camera System NS 4K resolution; 3840x2160

- 01. 4K Endoscope camera System NS 4K resolution; 3840x2160
- 02. Coupler, Water Proof
- 03. Monitor, 32" 4K
- 04. Trolley



- 01. 668 HD medical camera, 1080P
- 02. Coupler, Water Proof
- 03. 27" HD Monitor
- 04. Trolley, with multifunctional, shelf



- 01. 618 CCD camera
- 02. Coupler, F23, Water Proof
- 03. 21" Sony or 15" JVC Monitor
- 04. Pushcard, with multifunctional shelf
- 05. Monitor, 19" LCD
- 06. Printer, high resolution print
- 07. Food pedal
- 08. Video capture card, Philips chip highest quality of image, HD digital DVD level, dynamic magnification, dynamic filter processing
- 09. Software, can make DVD replay/record, print, many template/word/stock/photo stock, powerfull function of case management
- 10. Video teaching of the system, comprehensive teaching of the whole system
- 11. Automatic system recovery to factory setting in 20mins





Ginekologia

Endoscopy camera system

- 01. Full HD camera, LCD screen, 1080P, immersible camera head for disinfection
- 02. Water Proof coupler (Zoom coupler is optional)
- 03. 27 inch medical monitor
- 04. LED cold light source
- 05. Trolley with drawer and patch board



- 01. HD camera, 1080P, immersible camera head for disinfection
- 02. F14 Water proof coupler
- 03. 24inch medical monitor
- 04. LED cold light source
- 05. Trolley with drawer and patch board



Surveillance Monitor



- A 2502.37 Sony, 32 inch
- A 2502.38 Sony, 55 inch



- A 2502.24 24 inch, High definition, made in China
- A 2502.31 27 inch, High definition, made in China
- A 2502.39 32 inch, 4K



Ginekologia




Waginoskopy

Zestaw trzech wzierników z głowicami światłowodowymi, obturatorami w trzech rozmiarach i lupą powiększającą.



Kleszcze biopsyjne sztywne F1110.1
Kleszcze biopsyjne giętkie F1110



Amnioskopy

Nowa konstrukcja wziernika amnioskopowego, oświetlenie doprowadzone pod sam koniec wziernika. Tubus zintegrowany z głowicą światłowodową w trzech rozmiarach.



Wziernik amnioskopowy Ø 12 mm L= 200 mm



Wziernik amnioskopowy Ø 16 mm L= 200 mm



Wziernik amnioskopowy Ø 20 mm L= 200 mm

Szczypy biopsyjne do amnioskopu





Ginekologia

Light sources



M-OSW-LED-120-HCM-9887
Light source HCM 9887 LED 120W



M-OSW-LED-50-HCM-9671
Light source HCM 9671 LED 50W




M-OSW-HAL-150-1-XD 301-S
Light source XD 301, 150W, jednokanałowe



M-OSW-HAL-150-2-XD 301-S
Light source XD 301, 2 x 150W, STORZ



M-OSW-HAL-250-2-XD 301-S
Light source XD 301, 2 x 250W, STORZ



M-OSW-R-LED-3W-1
Portable LED light 3W

Endoscopic cameras



M-KAM-HD-UC-100
Endoscope camera UC-100 Full HD USB 3.0



M-KAM-CYMO-HD 1156
Endoscope camera CYMO HD 1156



M-KAM-CYMO-HD 6180
Endoscope camera CYMO HD 6180

Medical optica fiber

Adapter do kabla uniwersalnego firmy Optec, nakręcany od strony przyrządu.



Stryker



Wolf



Storz/Olympus

Adapter ksenonowy/halogenowy do kabla uniwersalnego firmy Optec, nakręcany od strony źródła światła.



Storz



Wolf



Olympus

M-SM-1800-35 -U-U
M-SM-2300-35 -U-U
M-SM-3000-35 -U-U
M-SM-3500-35 -U-U
M-SM-1800-48 -U-U
M-SM-2300-48 -U-U
M-SM-2500-48 -U-U
M-SM-3000-48 -U-U
M-SM-3500-48 -U-U

Medical opticalfiber Optec 1,8m $\varnothing$ 3,5Universal/Universal
Medical opticalfiber Optec 2,3m $\varnothing$ 3,5Universal/Universal
Medical opticalfiber Optec 3,0m $\varnothing$ 3,5Universal/Universal
Medical opticalfiber Optec 3,5m $\varnothing$ 3,5Universal/Universal
Medical opticalfiber Optec 1,8m $\varnothing$ 4,8Universal/Universal
Medical opticalfiber Optec 2,3m $\varnothing$ 4,8Universal/Universal
Medical opticalfiber Optec 2,5m $\varnothing$ 4,8Universal/Universal
Medical opticalfiber Optec 3,0m $\varnothing$ 4,8Universal/Universal
Medical opticalfiber Optec 3,5m $\varnothing$ 4,8Universal/Universal

$\varnothing$ 5.0< kabel nietypowy; wycena indywidualna



Ginekologia



Medical instruments sterylizing box




A0000.1
For Cystoscopy set 530x230x90mm



A0000.2
For Resectoscope set 420x275x90mm




A0000.4
For endoscope, autoclavable 390x71x36mm




A0000.5
Sterylizing box for Ureteroscopy
720x120x70mm




A0000.7
For Cystoscopy set 390x80x52mm nets



A0000.8
Sterylizing box for Ureteroscopy 660x140x90mm



A0000.9
300x140x90mm




A0000.10
500x110x120mm




F0000
For Electric hystera-cutter 560x275x140mm




F0000.1
Sterylizing box for Uterine manipulator 500x180x70mm nets



F0000.3
540x250x70mm



F0000.4
540x250x70mm



F0000.5
400x140x70mm



F0000.6A
400x110x45mm



F0000.8A
530x150x45mm



E0000A
For Laparoscopy 450x90x45mm



E0000.3
Sterylizing box for Laparoscopy 460x80x52mm nets




E0000.4
Reticular sterilizing box 100x100x55mm nets



E0000.11
550x250x100mm



190001000a
Endoscope Instruments in common use,
550x245x70mm



19.0001.000D
550x250x90mm




Ginekologia

Maintenance services

Usługi serwisowe:

- naprawa endoskopów
- naprawa światłowodów



Usługi serwisowe

WYSOKA JAKOŚĆ USŁUG - GWARANCJA - ATRAKCYJNE CENY - SZYBKI SERWIS

Firma **OPTEC** od wielu lat prowadzi serwis gwarancyjny i pogwarancyjny sprzętu medycznego wykorzystującego włókna światłowodowe. Nasze naprawy przeprowadzane są przez wykwalifikowanych pracowników przy użyciu specjalistycznego sprzętu i urządzeń. Dysponujemy ogromnym doświadczeniem, wiedzą oraz zapleczem technicznym, które pozwala na dokonanie profesjonalnych napraw.

Prowadzimy serwis i przeglądy techniczne następujących sprzętów:

- optyk endoskopowych sztywnych (artroskopy, laryngoskopy, otoskopy, sinuscopy, etc.),
- endoskopów giętkich (fiberoskopy, gastroskopy, kolonoskopy, nasofaryngoskopy, etc.),
- światłowodów medycznych oświetleniowych,
- światłowodów do laserów,
- elementów światłowodowych haków chirurgicznych, rektoskopów, amnioskopów, anoskopów, waginioskopów.

Dzięki naszemu zaangażowaniu możemy pochwalić się szeroką rzeszą stałych klientów.

To co nas wyróżnia na tle podobnych firm, to:

- bezpłatna wycena kosztów naprawy,
- stosowanie materiałów najwyższej jakości,
- krótkie terminy realizacji,
- gwarancja na wykonane czynności i zastosowane elementy.

**Jeżeli macie państwo jakieś pytania, chcecie się dowiedzieć czy możemy naprawić państwa sprzęt ?
Zapraszamy serdecznie do kontaktu telefonicznego lub mailowego z naszą firmą.**





Ginekologia




Maintenance services

Naprawa światłowodów


Specjalizujemy się w naprawach wszystkich typów światłowodów, w tym kabli do laserów i aparatury optoelektronicznej.

Włókna światłowodowe są najbardziej narażonym na uszkodzenia elementem aparatury, a ich zniszczenie uniemożliwia pracę całego systemu. Dlatego właśnie dbamy aby naprawy wykonywane były możliwie najszybciej.

Na przykład bardzo częsta usterka jaką jest uszkodzenie włókna przy konektorze naprawiana jest nie dłużej niż w dwa dni od dostarczenia.



Czolo światłowodu przed naprawą



Czolo światłowodu po naprawie



Wykonujemy regenerację osłony ochronnej lub jej wymianę na nową.

Istnieje możliwość instalacji nowego włókna światłowodowego.

Naprawiamy zniszczone końcówki bądź zmieniamy na inne.

Bezpłatnie określamy koszty naprawy

Stosujemy materiały i kleje najwyższej jakości

Na naprawiony sprzęt udzielamy 6 miesięcy gwarancji

OPTEC s.c. jest firmą z ponad 20-letnim doświadczeniem działającą w branży medycznej. Specjalizujemy się we wszelkiego rodzaju kablach światłowodowych do medycyny oraz optykach endoskopowych. Oferujemy sprzęt medyczny do artroskopii, ginekologii, laryngologii, laparoskopii, proktologii, urologii w bardzo konkurencyjnych cenach.



W naszej ofercie znajdują się: wszelkiego rodzaju światłowody do medycyny, źródła światła, oświetlacze czołowe, optyki endoskopowe sztywne i giętkie, wzierniki, tory wizyjne, kamery, unity, narzędzia chirurgiczne, haki operacyjne, lupy zabiegowe, wózki endoskopowe, pojemniki sterylizacyjne oraz wiele innych produktów. posiadamy prawie wszystkie rodzaje optyk sztywnych.

Dużym atutem naszej firmy jest szybki i atrakcyjny cenowo serwis wszystkich typów światłowodów, endoskopów sztywnych i giętkich (Storz, Wolf, Olympus, Stryker, Fujinon, Pentax, ACMI i inne), a także światłowodów specjalnych np. do lasera(lipoliza, metoda EVLT, urologia).

Optec s.c.

ul. Alberta 6E, 16-001 Księżyno k. Białegostoku

tel. 85 743 74 17, 85 878 16 10

e-mail: [medycyna @optec.pl](mailto:medycyna@optec.pl)

